

beyond[®] II Ultra Whitening Accelerator

Patented Ultrasonic Technology

Allows for the effective treatment of **two patients** at the same time.

Simultaneous*



Simultaneous



Ultrasound



Purifies



Adjustable



Digital



Ultra – the ultrasound function utilizes a cavitation process that permits the whitening gel to target the enamel more comprehensively and in a shorter time frame.

Simultaneous* - Allows for the effective treatment of two patients at the same time.*

Unique – Chic and modern design that will be easily embraced by any dental practice. Advanced design allows for a portable system, allowing for ease of mobility through different rooms within the practice.

Adjustable – Whitening option can be set to gentle, normal, or intense light output settings for a tailored treatment that best suits each individual patient.

Digital – Simple and quick patient preparation and system operation with digital display control board.

Purifies – Using negative ions, the built in air purifier reduces airborne pollutants that may exist within the practice, actively contributing to a cleaner and healthier environment.



The BEYOND® II Ultra whitening accelerator is a more portable option for in-office teeth whitening. The newest incarnation still allows for the simultaneous treatment of two patients. The tailored beam of blue light delivered by LEDs is only enhanced by the ultrasonic power that assists in providing superior whitening results in a short treatment time. Variable settings ensure that each patient receives a comfortable, positive whitening experience specifically catered to them. The three output settings of high, medium, or low on the digital display makes this specialized treatment possible. The lightweight portable design enables travel and ease of use from room to room. Utilizing negative ions, the built-in air purifier reduces airborne pollutants, actively contributing to a cleaner healthier environment.

Various colors available.

* Second LED sold separately.



beyond® International Inc.
— Global leaders in aesthetic dentistry. Products of USA
www.beyonddent.com

BEYOND USA Headquarters
711 Julie Rivers Dr.
Sugar Land, TX 77478 USA
TEL +1-281-277-4352
E-mail: info@beyonddent.com

BEYOND European Headquarters
Piotra Skargi 12
65-416 Zielona Gora, Poland
TEL: +48 512 401057
Email: europa@beyonddent.com

BEYOND Asia Pacific Headquarters
Chaoyang District, Beijing
China 100029
E-mail: asia@beyonddent.com

Questions? Call +1-855-956-6668 M-F 9:00 a.m. – 5:30 p.m. CST

YS00XCA0027,Rev.B