




BEYOND DENTAL PRODUCTS

SHOWCASE



60,700+
Treatment Locations
Worldwide



501,600+
Satisfied Dentists



Over 5 Million
Optimized Smiles



Since our inception in Sunnyvale California, BEYOND has focused on bringing comfortable, long lasting whitening results to smiles across the globe. Over 60,000 systems worldwide attest to our impact since 2003.

GLOBAL LEADER

of In-Office Teeth Whitening!

POLUS[®]
Advanced Ultra
Whitening Accelerator

Top Whitening System Award

15th Year in a row!



BEYOND patented ultrasonic technology cavitates the whitening gel to more comprehensively target the enamel. This allows us to achieve excellent whitening results under a shorter time frame.




**Touch User Interface (TUI)
Control Panel**



VOICE INSTRUCTIONS
**Remote Controlled Voice
Instruction Function**



**LED Curing Light /
Spot Whitening Device**




Using Ultrasonic Technology with our Teeth Whitening?
We like the sound of that!



POLUS[®]
Advanced Ultra
Whitening Accelerator

BEYOND® Polus® delivers the **Safest and Most Effective** whitening experience with customizable patient settings and a premium, ergonomic design for optimal performance.



Every BEYOND® Whitening Accelerator Comes With
A Complete Marketing Package

beyond® II Ultra
Whitening Accelerator

Whitening Made Simple and Delicate!



- High-Intensity LED with Ultrasound Technology
- Portable & Compact Design
- Super-Efficient Treatment Time
- Affordable Pricing

* Second LED sold separately.



beyond[®] II Ultra
Whitening Accelerator



ULTRASOUND



ADJUSTABLE



DIGITAL



PURIFIES

Patented Ultrasonic Technology
Allows for the effective treatment of **two patients**
at the same time.



Professional in-office Whitening Treatment

Every BEYOND® II Whitening Accelerator Comes With
A Complete Marketing Package



BEFORE AFTER



BEYOND InPulse

The “**Nd:Yag**” you are looking for!

INPULSE


INPULSE DENTAL LASER
(Nd:YAG) InPulse Touch 100





INPULSE

INPULSE DENTAL LASER (Nd:YAG) InPulse Touch 100



Experience precision and comfort with BEYOND InPulse Nd:YAG laser technology. With its ability to penetrate tissue without being easily absorbed by water, and the use of pulse oscillation for precise heat penetration control, treatments can be performed with less impact and risk. In most cases, no anesthesia is required for a pain-free experience.



BEYOND® Series
Teeth Whitening Gel
is Recommended by **Dental Advisor** as

The Best



In-Office Teeth Whitening Kits





BY-MX105-US

FORMULATED TO WORK WITH ANY WHITENING

BEYOND® MAX5® Five-Patient Treatment Kit

BEYOND's Max5 is our most popular In-office teeth whitening kit. This 35% hydrogen peroxide (H_2O_2) kit includes the essentials for chairside teeth bleaching of five patients. A typical treatment consists of three consecutive cycles of 10 minutes each, in conjunction with light acceleration. The dual barrel 35% syringe in this kit includes gel for three cycles.

FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.

Package Contents

- 5 Max5 gel dual-barrel syringes (2.6g / 2.2 mL)
- 5 BlueSeal single syringes (2g / 1.6mL)
- 5 Mixing tips
- 5 Plastic tips
- 5 Metal tips
- 5 Face protection cloth
- 1 Instruction manual
- 5 Vita Shade Guides
- 5 After-treatment guides



BY-SL101

OUR MOST ECONOMICAL, SINGLE PATIENT WHITENING KIT

BEYOND® SOLO® Single Treatment Kit

This 35% hydrogen peroxide (H_2O_2) In-office teeth whitening kit includes the essentials for chairside teeth bleaching of one patient. A typical treatment is three consecutive cycles of 10 minutes each, combined with light acceleration. The dual barrel 35% syringe in this kit includes enough gel for three cycles. This professional in-office kit is the single patient version of our Max5 kit.

Package Contents

- 1 Solo gel dual-barrel syringe (2.6g / 2.2 mL)
- 1 BlueSeal single syringe (2g / 1.6mL)
- 1 Mixing tip
- 1 Plastic tip
- 1 Metal tip
- 1 Face protection cloth
- 1 Instruction manual
- 1 Vita Shade Guide
- 1 After-treatment guide



BY-PD202

SINGLE PATIENT, COMPREHENSIVE WHITENING PACKAGE

This 35% hydrogen peroxide (H_2O_2) In-office teeth whitening kit includes everything needed for one patient. The COMPLETE professional dental grade kit is a comprehensive choice as it combines both chairside and take-home products in one package. Designed to treat most common teeth bleaching cases, a typical treatment is three consecutive cycles of 10 minutes each combined with light acceleration.

BEYOND® COMPLETE® Professional Teeth Whitening Kit with Tray Whitening Material

The dual barrel 35% syringe includes enough gel for four cycles. The two CoreWhite syringes contain a gentle 6% H_2O_2 formula that you can send home with your patient for at home touch ups, as well as material for creating tailored trays. In conjunction with the Polus Advanced system, the Dental Advisor awarded this kit as the 2016 recipient of the best Whitening: In-Office.

Package Contents

- 1 Beyond single dual-barrel syringe 35% (4g / 3.5mL)
- 2 CoreWhite home gel syringes (mint) 6% (2.36g / 2.2mL)
- 1 BlueSeal single syringe (2g / 1.6mL)
- 1 Mixing tip, 1 Plastic tip, 1 Metal tip
- 1 Suction tip
- 1 Cheek retractor,
- 1 Face protection cloth
- 1 Jar of lip cream (2.4g / 2mL)
- 2 Sheets EVA impression tray material
- 1 Instruction manual
- 1 Vita Shade Guide
- 1 After-treatment guide



BY-PD205

OUR MOST ECONOMICAL, FULL-SIZED WHITENING KIT

This 35% Hydrogen Peroxide (H_2O_2), In-office teeth whitening kit includes material for five patients. The Five-Patient Economy kit is a professional dental grade package designed to treat most common teeth bleaching cases. A typical treatment is three consecutive cycles of 10 minutes each, combined with light

BEYOND® Five-Patient Economy Kit

acceleration. The dual barrel 35% syringe in this kit includes enough gel for four chairside cycles.

Package Contents

- 5 Beyond single dual-barrel syringes (4g / 3.5 mL)
- 5 BlueSeal single syringe (2g / 1.6mL)
- 5 Mixing tips
- 5 Plastic tips
- 5 Metal tips
- 5 Cheek retractors
- 5 Jars of lip cream (2.4g / 2mL)
- 5 Face protection cloth
- 1 Instruction manual
- 5 Vita Shade guides
- 5 After-treatment guides



BY-OS201

ADVANCED 12% H₂O₂ WHITENING KIT WITH LASTING PEN

Introducing the upgraded OSMO Plus Five-Patient Whitening Kit with 12% hydrogen peroxide gels for faster, more effective results, allowing patients to achieve brighter smiles in just 36 minutes. Our advanced, safe-to-use formula requires no gingival protection, simplifying the whitening process. The kit, sufficient

BEYOND® OSMO® Plus Five-Patient Professional Kit

for five patients, includes a dual barrel syringe for chairside treatments and premium whitening pens for home maintenance, ensuring long-lasting, dazzling smiles.

Package Contents

- 5 Osmo Plus gel dual-barrel syringes (2.6g / 2.2ml)
- 5 Osmo teeth whitening pens (3.2g each)
- 5 Mixing tips
- 5 Face protection cloth
- 1 Instruction manual
- 5 Vita Shade Guides
- 5 After-treatment guides



BY-OS101

6% H₂O₂ PREMIUM WHITENING FORMULA WITH POST-TREATMENT WHITENING PEN

This 6% Hydrogen Peroxide (H₂O₂) Teeth Whitening Kit is designed to deliver professional results in just 36 minutes. The OSMO Five -Patient Whitening Kit features a safe, advanced formula gel that can be used in-office without gingival protection. This kit includes material for five patients. A typical treatment is three consecutive cycles of 12 minutes each, designed to work in conjunction with light acceleration. The dual barrel 6% syringe in this kit includes enough teeth bleaching gel for

BEYOND® OSMO® Five-Patient Professional Kit

three chairside cycles. The Osmo Kit also includes five premium formula teeth whitening pens that patients can take home and start using immediately to keep stains and darkening at bay while adding a rich, glossy tooth finish.

Package Contents

- 5 Osmo gel dual-barrel syringes (2.6g / 2.2ml)
- 5 Osmo teeth whitening pens (3.2g each)
- 5 Mixing tips
- 5 Face protection cloth
- 1 Instruction manual
- 5 Vita Shade Guides
- 5 After-treatment guides



BY-OS203 / BY-OS203A

SALON-GRADE WHITENING
PERFECTION IN ONE
COMPREHENSIVE PACKAGE.

BEYOND® WhiteSpa®
Single-Patient Kit

Enhance your salon with the WhiteSpa Whitening Kit, featuring a unique 6% /12% H₂O₂ gel in a dual barrel syringe and designed to work with any BEYOND whitening accelerator for superior results. This innovative system offers a cutting-edge approach to teeth whitening, ensuring remarkable effectiveness and client satisfaction with a dazzling smile.

Package Contents

- 1 Dual-barrel syringe of WhiteSpa premium formula whitening gel (2.6g / 2.2mL)
- 1 Mixing tip
- 1 Osmo teeth whitening pen (2.3g / 1.9mL)
- 1 Bib
- 3 Vitamin E cotton swabs
- 4 Finger wipes
- 1 Face protection cloth
- 1 VITA Shade Guide
- 1 Instruction manual
- 1 After-treatment guide



BY-BW101US

SPECIAL DESIGNED FORMULA
FOR TETRACYCLINE AND
SEVERE STAINS

BEYOND® Blast®
Pre-Treatment Extra-Strength
Whitening Paste

The BEYOND BLAST is a professional dental grade pretreatment paste designed to effectively jump-start the teeth whitening process. BEYOND BLAST's 35% hydrogen peroxide (H₂O₂) formula is especially effective as a pre-treatment for tetracycline and age-related stains. The process is simple: after patient preparation, mix and apply paste, and introduce light acceleration source.

Package Contents

- 1 Bottle of 35% H₂O₂ (2mL)
- 1 Jar of catalyst
- 1 Spatula
- 1 Application brush
- 1 Face protection cloth
- 1 Instruction manual
- 1 After-treatment guide



REBUILD —hydroxyapatite

STRENGTHEN —potassium nitrate

WHITEN —hydrogen peroxide, calcium peroxide

The BEYOND® Solution range of kits work to deliver excellent results and with added desensitizing gel for those patients more susceptible to sensitivity issues. Our BEYOND® RS gel is designed to restore and desensitize teeth after treatment, ensuring comfort and longevity of treatment results.


Our BEYOND® Solution series of whitening kits are revolutionary, combining high percentage in-office whitening kits with a specially formulated RS gel to reduce known sensitivity of some patients. This kit includes the same original dual-barrel syringe with 35% hydrogen peroxide (H_2O_2) that our partners are familiar with and our specially formulated RS Gel 2 that contains desensitizing ingredients. This combination of potassium nitrate and hydroxyapatite offers enamel regeneration and restoration. When combined, the strength of both gels reduce the hydrogen peroxide concentration to 6% and helps aid in the reduction of sensitivity as well as remineralization of the enamel.

Using the RS gel alone enables patients to overcome any hypersensitivity in a very short time after treatment and does not require the use of the accelerator. The end result is excellent whitening coupled with little to no sensitivity for those undergoing the treatment.

• BlueSeal Gingival Protection sold separately



With our EU headquarters,
we ensure fast delivery across
Europe within days.



BEYOND® Solution

Range of Kits

The Solution Kits Series is designed for markets like Europe, complying with the 6% H_2O_2 limit while ensuring safe and effective whitening.





BY-MS105

REBUILD, STRENGTHEN, AND
WHITEN FORMULA

BEYOND® MAX5® Solution Treatment Kit

Max Solution delivers the essentials for chairside teeth bleaching of five patients. The dual barrel 35% W Gel syringe for this kit includes enough gel for three cycles. RS gel can be then applied after treatment. BEYOND® RS gel is designed to restore and desensitize the teeth after treatment, ensuring comfort and longevity of treatment results.

Package Contents

- 5 Dual-barrel syringes of Max5 Solution W gel (2.6g / 2.2mL)
- 5 Mixing tips of Max5 Solution W gel syringe
- 5 Pouches of Max5 Solution RS gel (11g / 9.5mL)
- 5 Spatulas, 5 Brushes, 5 Mixing jars
- 5 Face protection cloths
- 5 Vita shade guides
- 5 After-treatment guides
- 1 Instruction manual



BY-SS101

OUR MOST ECONOMICAL,
SINGLE PATIENT WHITENING KIT

BEYOND® SOLO® Solution Treatment Kit

BEYOND® Solo® Solution is a single patient version of the Max5 Solution Kit. Solo® Solution consists of a single W Gel dual barrel syringe of 35% H₂O₂, which is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredients such as hydroxyapatite.

Package Contents

- 1 Dual-barrel syringe of Solo Solution W gel (2.6g / 2.2mL)
- 1 Mixing tip for Solo Solution W gel syringe
- 1 Pouch of Solo Solution RS gel (11g / 9.5mL)
- 1 Face protection cloth
- 1 Spatula
- 1 Brush
- 1 Vita shade guide
- 1 Instruction manual
- 1 After-treatment guide



BY-MS102

ECONOMICAL FULL-SIZE, TWIN-PACK WHITENING KIT

BEYOND® MAGNA™ Solution Treatment Kit

BEYOND® Magna Solution is a full size two patient treatment kit, a simplified and enlarged version of the Max5 Solution kit. Each W Gel syringe contains enough gel for initial treatment as well as additional cycles as necessary or for a standalone touch up treatment .

Staff can alternatively use the provided spatula and Jar to apply RS Gel directly to teeth for special sensitivity cases.

Package Contents

- 2 Magna Solution W gel 1 dual-barrel syringes of 35% H_2O_2 (4.0g / 3.5 mL)
- 2 Mixing tips of Magna Solution W gel dual-barrel syringe
- 2 Pouches of Magna Solution RS gel (16g / 13.8mL)
- 2 Brushes, 2 Spatulas
- 2 Face protection cloths
- 2 Vita shade guides
- 1 Instruction manual



BY-BL101

SPECIAL DESIGNED FORMULA FOR TETRACYCLINE AND SEVERE STAINS

BEYOND® Blast® Solution Treatment Kit

BEYOND® Blast Solution is a professional dental grade pretreatment paste that effectively jump-starts the teeth whitening process. Especially effective as pre-treatment for tetracycline and age-related stains. This treatment can be followed with any other BEYOND In-office product.

Package Contents

- 1 Jar of Blast® Solution silicon dioxide catalyst powder (0.7g)
- 1 Bottle of Blast® Solution 35% hydrogen peroxide aqueous solution (2mL)
- 1 Pouch of Blast® Solution RS gel (16g / 13.8mL)
- 1 Spatula, 1 Brush
- 1 Face protection cloth
- 1 Vita shade guide
- 1 Instruction manual

BEYOND® CoreWhite® Take-Home Whitening Kit

- 30 minutes a day for 14 days – no messy overnight whitening
- Visibly whiter teeth in five days
- Gentle hydrogen peroxide formula minimizes risk of sensitivity (6% H₂O₂ gel)
- EVA material included for a comfortable customized tray / A set of universal tray (two options)
- Mint and Orange flavors
- CoreWhite® Refill Kits also available



Our gentle 6% hydrogen peroxide (H₂O₂) formula promises professional results from the comfort of your home. The CoreWhite take-home tray teeth whitening system offers a quick and inexpensive treatment lasting only 30 minutes a day while not requiring custom trays, providing convenience for a beautiful white smile.

Package Contents

- 6 syringes CoreWhite gel (2.36g / 2.2 mL)
- 2 Sheets of EVA impression material
- 1 CoreWhite tray box
- 1 Manual booklet
- 1 Carrying case

Take - Home Tray Whitening
Available through the Dental Office



CoreWhite® Refill Kits also available



BEYOND® OSMO®

Teeth Whitening Pen Twin Pack

This pack features two teeth whitening pens, each with enough hydrogen peroxide (H_2O_2) gel to last for 30 days. Use the pen for just minutes daily as maintenance following a professional teeth whitening service or as a stand alone at-home teeth whitening product.



osmo®

Package Contents

2 Osmo® whitening pens (3.2g / 2.8 mL)


- Provides a renewed white smile every day
- Safe for daily use and whitening maintenance
- Gentle brush-on pen for easy application
- Can be used alone or in conjunction with all
- BEYOND® whitening procedures
- Enough material for 30 days
- Mint flavored 6% H_2O_2 gel for a white and fresh smile

Revolutionize your teeth in just 1 Hour

BEYOND® Gemini® Teeth Whitening Accelerator Home Edition



gemini®



- Next generation LED for maximum acceleration
- No sensitivity
- Removes years of stains caused by coffee, red wine and smoking
- 6 to 8 shade improvements
- Embedded magnets line light to tray for a perfect fit every time
- Simply rechargeable via USB or outlet
- Clinically proven to kill germs
- Refill kits available for continued maintenance

Package Contents

- 1 LED light unit
- 5 Syringes of Gemini® whitening gel (1.9g / 1.7mL)
- 5 Tips of Gemini® whitening gel syringe
- 1 Gemini® tray
- 1 power unit with rechargeable lithium battery
- 1 charge cord
- 1 USB / electrical adaptor
- 1 Bleaching shade guide
- 1 Instruction manual

PEARL WHITE®

Whitening Toothpaste



Green Tea Formula

- Includes camelia sinensis extract (green tea extract) that:
 - Act as an antioxidant. Great for the health of gum tissue
 - Inhibit odor and lasting fresh breath
- Unique pearl powder formula leaves teeth shiny and sparkling white



Sensitivity Formula

- Soothes irritated nerves in your teeth
- Protects against the causes of hypersensitivity
- Active ingredients : Sodium Fluoride and Potassium Nitrate
- FDA registered



Advanced Formula

- Long-lasting, minty-fresh taste
- Unique pearl powder formula leaves teeth shiny and sparkling white



Advanced Daily Whitening Maintenance

beyond[®]

Nature Series Toothpaste

*Think
Think*

*Fresh,
Nature.*



Fresh & Healthy

Daily Whitening Maintenance

BEYOND® is proud to announce our newest additions to our award-winning line of whitening toothpaste. In addition to our legacy Pearl White whitening toothpaste, we are now pleased to introduce the next-generation formulas.



beyond® ROSEMINT TOOTHPASTE

- Gently Whitens
- Soothes Gums
- Moisturizes Mouth
- Freshens Breath



beyond® MATCHA TOOTHPASTE

- Gently Whitens
- Prevents Plaque Buildup
- Freshens Breath
- Improves Gum Health



beyond® Activated Coconut CHARCOAL TOOTHPASTE & SEA SALT

- Deeply Whitens
- Fights Cavities
- Reduces Plaque
- Soothes Gums



BEYOND® Premium Bamboo Toothbrush

The BEYOND Premium Bamboo toothbrush is **eco-friendly** with a fully biodegradable bamboo handle. The soft vegan-friendly, **BPA-free** nylon bristles allow you to clean your teeth more thoroughly.



Package Contents

6 BEYOND® bamboo toothbrushes



- **ECO-FRIENDLY** - If you're looking for an environmentally friendly toothbrush option, our bamboo toothbrush is the perfect choice. It's made of sustainable bamboo, and the quality is just like a regular plastic toothbrush. Our packaging is also 100% recyclable.
- **BPA-FREE SOFT BRISTLES** - Soft bristles that clean teeth without hurting gums, and are crafted to never shed.
- **ERGONOMICALLY DESIGNED HANDLE** - fits comfortably in your hand, easy to grip and control.
- **PERFECT FOR FAMILIES** - Our package of 6 toothbrushes with 3 different colors enable easily to tell every toothbrush apart, very suitable for family members. Each toothbrush lasts up to 6 months of use, making this a cost-effective way to keep the family brushing.

High Quality Dental Accessories

Your Whitening

ONE-STOP SHOP



Convenient, Gingival Protection Syringes

BEYOND® BlueSeal® Gingival Protection

Package Contents

2 / 5 / 10 BlueSeal Syringes (2.36g / 2.2 mL)

2 / 5 / 10 Metal tips

2 / 5 / 10 Plastic tips

BlueSeal™

BlueSeal creates a precise gingival barrier required during the use of dental grade teeth whitening gels. BEYOND's gingival protection comes with 2 different types of tips and cures in just minutes. For use with any brand of teeth bleaching.

- Gingival protection barrier for dental use
- Easy-to-use syringe for precision application
- Cures within seconds and stays in place
- Simple post-treatment removal



FlexiBlue®

BEYOND® FlexiBlue® Comfort Cheek Retractors

- Easy to insert, available in two size
- Convenient tongue guard
- No sharp edges

BEYOND's custom designed Flexiblu cheek retractors are designed to work with any accelerator. The comfortable design holds lips away from enamel for the duration of teeth whitening processes. The innovative design includes a built in tongue depressor. Available in two sizes: medium or large.



Medical Grade, Disposable Cheek Retractors



BEYOND® FlexiPro™,
More Comfortable than Ever!

FlexiPro™

The BEYOND® FlexiPro™ Cheek Retractor is an in-office dental device intended for increasing the work area of the mouth during teeth whitening procedures.

The BEYOND® FlexiPro™ Cheek Retractor contains a pair of channeled edges that can be comfortably fitted around the lips of the patient and a bridge connecting the edges of the device. The bridge features a tongue guard that holds back the patient's tongue during procedures, while doubling as a bite cushion for the patient to rest their teeth on.

BEYOND® FlexiPro™ Cheek Retractor is designed to enlarge the vestibular fundus and increase the work area of the anterior 16-20 teeth, as well as provide protection to the lip area. The softer material and elasticity provides patients with a safe and comfortable treatment condition. The device can be cold sterilized.

BEYOND® FlexiPro™ Cheek Retractor





BEYOND® Protective Eye Goggles

BEYOND® Protective Eye Goggles safely and reliably absorb all blue light as protection from curing lights, computer screens and strong light/ultraviolet radiation. The goggles are also strong protection for eyes and skin from liquid they may be typically used in dental operations. These goggles are ergonomic to easily fit the face with no metal parts, making the user more comfortable.

- Metal frame
- Polycarbonate lens
- Anti-Fog Coating coating
- Made with 100% reinforced polycarbonate
- Complies with American standard ANSI Z87.1 level 1 impact resistance standard and European standard CE EN166 test standard
- Wing and eyebrow guard design block particles, dust and/or liquid originating from sides
- Orange lens, blue light absorption, Medical Dental Protective Safety Glasses for Dentist, Hygienist, Assistant, Patients, Lab personnel
- Adjustable-length temples, molded-in nose bridge



Make your patient and you
A Comfort and Safety treatment experience!

beyond®

Quick, Secure And precise shade determination

BEYOND® 3D

Teeth Whitening Shade Guide

- 20 shades with mirror case, more professional. Laser engraved numbers are easy to read and longer lasting use
- The shade tabs can be heat-sterilized in the autoclave up to a max temperature of 140 or cold-sterilized by using a suitable disinfectant
- The shade guide base is only suitable for surface disinfection
- Note: the plastic parts of the shade guide are not resistant to monomers



beyond[®] *International Inc.*
beyonddent.com

BEYOND USA Headquarters

711 Julie Rivers Dr.
Sugar Land, TX 77478 USA
TEL +1-281-277-4352
Toll Free +1-855-956-6668
E-mail: info@beyonddent.com

BEYOND European Headquarters

Plotra Skargl 12
65-416 Zielona Gora, Poland
TEL +48 512 401057
E-mail: sales@beyond-europe.biz

BEYOND Asia Pacific Headquarters

Chaoyang District, Beijing, China 100029
E-mail: asia@beyonddent.com