



### BEYOND® MAX5® Solution Treatment Kit

The BEYOND® Max5® Solution Treatment Kit pairs our popular in-office whitening kit with RS Gel to minimize sensitivity. It includes a W Gel 1 dual-barrel syringe with 35% hydrogen peroxide and RS Gel 2 with potassium nitrate and hydroxyapatite. Together, they reduce hydrogen peroxide to 6%, lower sensitivity, and promote enamel remineralization. Designed for five patients, each treatment involves three 15-20 minute chairside cycles with light acceleration. RS Gel alone reduces post-treatment sensitivity, providing excellent whitening results with minimal discomfort.

BY-MS105



### BEYOND® SOLO® Solution Treatment Kit

BEYOND® Solo® Solution was created with patients with known cases of sensitive teeth in mind. Solo® Solution consists of a W Gel of 35% H<sub>2</sub>O<sub>2</sub> dual barrel syringe, which is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredients such as hydroxyapatite. This reduces the concentration of hydrogen peroxide to below 6% and helps aid in reduction of sensitivity as well as remineralization of the enamel. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment.

BY-SS101



### BEYOND® MAGNA™ Solution Treatment Kit

BEYOND® Magna™ Solution was created with patients with known cases of sensitive teeth in mind. Magna™ Solution consists of a W Gel of 35% H<sub>2</sub>O<sub>2</sub> dual barrel syringe, which is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredients such as hydroxyapatite. This reduces the concentration of hydrogen peroxide to below 6% and helps aid in reduction of sensitivity as well as remineralization of the enamel. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment.

BY-MS102



### BEYOND® Blast® Solution Treatment Kit

BEYOND® Blast® Solution is designed for patients with severe or tetracycline stains. It combines a 35% H<sub>2</sub>O<sub>2</sub> W Gel with a silicon dioxide catalyst powder to form a wet paste, which is then mixed with RS Gel containing desensitizing potassium nitrate and restorative hydroxyapatite. This reduces the hydrogen peroxide concentration to below 6%, minimizing sensitivity and promoting enamel remineralization. The result is effective whitening with minimal discomfort.

BY-BL101



## **beyond®** Premium Whitening Solutions

# The Best In-Office Teeth Whitening Kits



BEYOND USA Headquarters  
711 Julie Rivers Dr.  
Sugar Land, TX 77478 USA  
TEL +1-281-277-4352  
E-mail: info@beyonddent.com

BEYOND European Headquarters  
Piotra Skargi 12  
65-416 Zielona Gora, Poland  
TEL: +48 512 401057  
Email: europe@beyonddent.com

BEYOND Asia Pacific Headquarters  
Chaoyang District, Beijing  
China 100029  
E-mail: asia@beyonddent.com

Questions? Call +1-855-956-6668 M-F 9:00 a.m. – 5:30 p.m. CST

YSOOXCA0026,Rev.B

www.beyonddent.com



# **beyond®** Premium Whitening Solutions



## **BEYOND® MAX5®** Five-Patient Treatment Kit

BEYOND's Max5 is our most popular In-office teeth whitening kit. This 35% hydrogen peroxide ( $H_2O_2$ ) kit includes the essentials for chairside teeth bleaching of five patients. A typical treatment consists of three consecutive cycles of 10 minutes each, in conjunction with light acceleration. The dual barrel 35% syringe in this kit includes gel for three cycles.

FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.

BY-MX105



## **BEYOND® SOLO®** Single Treatment Kit

This 35% hydrogen peroxide ( $H_2O_2$ ) In-office teeth whitening kit includes the essentials for chairside teeth bleaching of one patient. A typical treatment is three consecutive cycles of 10 minutes each, combined with light acceleration. The dual barrel 35% syringe in this kit includes enough gel for three cycles. This professional in-office kit is the single patient version of our Max5 kit.

BY-SL101



## **BEYOND® Simplex™** Five-Patient Treatment Kit

Unlock a dazzling smile with the Simplex Teeth Whitening Kit—a powerful 37% hydrogen peroxide formula that delivers stunning results in just 45 minutes, no light needed! Its user-friendly dual-barrel syringe with a dynamic mixing tip ensures effortless application, while the smooth gel texture allows precise placement and easy removal. Perfect for dental offices seeking professional whitening without extra equipment, Simplex offers a hassle-free, high-impact solution.

BY-SP106



## **BEYOND® COMPLETE®** Professional Teeth Whitening Kit with Tray Whitening Material

This 35% hydrogen peroxide ( $H_2O_2$ ) in-office teeth whitening kit includes everything needed for one patient. The COMPLETE professional-grade kit combines chairside and take-home products, designed to treat common teeth bleaching cases. A typical treatment involves three 10-minute cycles with light acceleration. The dual-barrel 35% syringe provides gel for four cycles. The two CoreWhite syringes, containing a gentle 6%  $H_2O_2$  formula, are suitable for at-home touch-ups and come with material for tailored trays.

BY-PD202



## **BEYOND® Five-Patient Economy Kit**

This 35% Hydrogen Peroxide ( $H_2O_2$ ), In-office teeth whitening kit includes material for five patients. The Five-Patient Economy kit is a professional dental grade package designed to treat most common teeth bleaching cases. A typical treatment is three consecutive cycles of 10 minutes each, combined with light acceleration. The dual barrel 35% syringe in this kit includes enough gel for four chairside cycles.

BY-PD205



## **BEYOND® OSMO®** Five-Patient Professional Kit

This 6% Hydrogen Peroxide ( $H_2O_2$ ) Teeth Whitening Kit delivers professional results in 36 minutes. The OSMO Five-Patient Whitening Kit features a safe, advanced gel formula suitable for in-office use without gingival protection. It includes enough material for five patients, with a typical treatment involving three consecutive 12-minute cycles using light acceleration. The 6% dual-barrel syringe provides gel for three chairside cycles. Additionally, the kit includes five whitening pens that patients can take home to prevent stains and maintain a glossy finish.

BY-OS101



## **BEYOND® OSMO® Plus** Five-Patient Professional Kit

Introducing the upgraded OSMO Plus Five-Patient Whitening Kit with 12% hydrogen peroxide gels for faster, more effective results, allowing patients to achieve brighter smiles in just 36 minutes. Our advanced, safe-to-use formula requires no gingival protection, simplifying the whitening process. The kit, sufficient for five patients, includes a dual barrel syringe for chairside treatments and premium whitening pens for home maintenance, ensuring long-lasting, dazzling smiles.

BY-OS201



## **BEYOND® Blast®** Pre-Treatment Extra-Strength Whitening Paste

The BEYOND BLAST is a professional dental grade pretreatment paste designed to effectively jump-start the teeth whitening process. BEYOND BLAST's 35% hydrogen peroxide ( $H_2O_2$ ) formula is especially effective as a pre-treatment for tetracycline and age-related stains. The process is simple: after patient preparation, mix and apply paste, and introduce light acceleration source.

BY-BW101

# Teeth Whitening One-Stop Shop