



Top Whitening System
Award **15th Year** in a row!




— Global leaders in aesthetic dentistry. Products of USA
www.beyonddent.com

BEYOND USA Headquarters
711 Julie Rivers Dr.
Sugar Land, TX 77478 USA
TEL +1-281-277-4352
E-mail: info@beyonddent.com

BEYOND European Headquarters
Piotra Skargi 12
65-416 Zielona Gora, Poland
TEL: +48 512 401057
Email: europe@beyonddent.com

BEYOND Asia Pacific Headquarters
Chaoyang District, Beijing
China 100029
E-mail: asia@beyonddent.com


Questions? Call +1-855-956-6668 M-F 9:00 a.m. – 5:30 p.m. CST

BEYOND patented ultrasonic technology
allows you to achieve excellent whitening
results under a shorter time frame.

POLUS[®]
Advanced Ultra
Whitening Accelerator

Using Ultrasonic Technology with our Teeth Whitening?

We like the sound of that!



LightBridge® Technology + Ultrasonic Technology



VOICE INSTRUCTIONS

Remote Controlled Voice Instruction Function



Touch User Interface (TUI) Control Panel

Beyond® Polus® Advanced Ultra Teeth Whitening System. the Fifteen Time Winner of the Top In-office Whitening Products (The Dental Advisor 2011-2025), BEYOND patented ultrasonic technology cavitates the whitening gel to more comprehensively target the enamel. This allows us to achieve excellent whitening results under a shorter time frame.

The LightBridge® is an innovation created to optimize the full whitening capabilities of halogen power while eliminating the need to block out the oral cavity. The LightBridge® consists of 30 specially coated lens that blocks all UV and harmful heat, coupled with over 200,000 optical fibers that deliver the full halogen output to both arches of the patient. This is the key to being able to deliver superior results in a shorter chairside time.




Unprecedented Innovation and Flexibility in Whitening Technology

- Upgraded light filtration system containing over 200,000 optical fibers, to filter out all harmful heat and ultraviolet (UV) light
- Dual-arch power whitening option can be set to High, Medium, or Low light output settings for a tailored treatment that best suits each individual patient
- Single-tooth whitening option included for precision, tailored treatment results
- Smart Remote Control that allows patient awareness of procedure time
- Built-in passive air ionizer emits negatively charged Built-in air purifier for a cleaner and healthier environment.

Innovative and Power-Packed System

- Improved motorized arm that allows for yet easier adjustment of the accelerator head for a more comfortable whitening procedure
- State of the art, full colored, touch screen with newly added multilingual feature
- Built-in instruction mode, that allows for a customized whitening procedure based on the client's needs
- Attached LED curing light, featuring a dual wavelength output to cure any dental composite on the market
- Radial arm with 360 range of motion
- Motorized height adjustment can be controlled from the TUI control panel and accelerator head
- Advanced temperature cooling and monitoring system designed to prevent system from overheating, enabling high volume whitening sessions without interruption.

TUI Control Panel



beyond®
POLUS®
Advanced Ultra